

Simple to use

Treatment for Hands



Treatment for Feet



Treatment for Underarms



The HIDREX therapy

Over 40 years ago, HIDREX developed with the largest university clinics in Germany local therapy against excessive sweating. In 1987, HIDREX became the first company to manufacture professional devices for the treatment of hyperhidrosis. The HIDREX therapy is now a scientifically recognized treatment method worldwide and is also known as tap water iontophoresis (TWI).

- Clinically tested with 98% success rate
- Fully automatic therapy procedure
- Automatic polarity change
- Powerful devices, up to 60V and 30mA
- ASE – anti-shock electronic
- Protection against over-treatment
- Immersion monitoring
- 100% German manufactured



Direct line to customer service!

Simply scan the QR Code and contact our customer service.

You will receive only answers to your questions and no newsletters or advertising.



HIDREX GmbH
Otto-Hahn-Str. 12
42579 Heiligenhaus
Germany
Fon: +49 (0) 2056 / 9811-0
Email: info@hidrex.de



HIDREX[®]

THERAPY

We stop sweating.

Treatment against excessive sweat!



Hands



Feet



Underarms



Face



www.hidrex-therapy.com

Hyperhidrosis

The term hyperhidrosis refers to excessive sweating that goes beyond the body's need for thermal regulation. People with hyperhidrosis sweat excessively, regardless of physical exertion or ambient temperature.

This dysfunction can lead to significant mental, social, and occupational limitations and can thus cause a marked reduction in the quality of life. About 5% of the population suffer from localized, including in the hands, feet and underarms. This disease is wrongly considered by many to be untreatable.



For more information visit
www.hidrex-therapy.com

Tap water iontophoresis (TWI)

Iontophoresis is a treatment method that has been scientifically recognized for many years and can result in normally dry skin without medication or surgery.

In comparison to many other therapeutic methods, iontophoresis has no negative long-term effects. The procedure is easy to perform, safe, inexpensive, and can be used even on teenager and children.

In addition, beyond the purchase of an iontophoresis device, there are no further follow-up costs, so that the treatment is very affordable of the long term.

Additional
accessories
available!



CE0633 
Cleared as a Class 2
Medical Device (510k)



HIDREX has distributors all over the world. Find your local supplier at hidrex-therapy.com/buy. For information on renting or buying the devices, contact your local supplier or contact us directly. Many insurance companies reimburse the cost of the device.



We stop sweating.

In HIDREX-therapy, a mild DC or pulse current is passed through the affected body zones. The flow of this electrical current is ensured by a water bath and leads to a sweat reduction on the areas being treated.

HIDREX treatment is administered by hand and foot baths or by treatment of the underarms, face, neck, or back with special sponge applicators.

During the initial treatment phase, you should, if possible, treat daily for 15-minutes. After about 10 treatments, your sweat levels should normalise. Afterwards, 1-3 treatments a week are usually sufficient to maintain the results. The therapeutic effect is reversible, so if treatments stop, the sweating will resume after some time. Therefore, the purchase of a home therapy device is recommended.